




T.T.A. Kit



TTA Box - Closed



TTA Box - Lid Removed



TTA Tray

Code:	Description:	Price:
<u>TTA Box</u>	For 2.4, 2.7 & 3.5 mm Screws, Instruments & Drill Bits - Empty 292 mm L X 145 mm W X 110 mm D	\$650.00

A case will hold the following:

Screws

- 2.4 mm - 10 each of 10, 12, 14, 16, 18, 20, 22, 24, 26, 28, 30, 32, 34, 36, 38, & 40 mm lengths (160 total)
- 2.7 mm - 10 each of 10, 12, 14, 16, 18, 20, 22, 24, 26, 28, 30, 32, & 34 mm lengths (130 total)
- 3.5 mm - 10 each of 10, 12, 14, 16, 18, 20, 22, 24, 26, 28, 30, 32, 34, 36, 38, & 40 mm lengths (160 total)

Additional Items

- Drill Bits - Jacobs Style - 1 each 1.8 mm, 2.0 mm, & 2.5 mm
- Drill Bits - Quick Coupling Style - 1 each 1.8 mm, 2.0 mm, & 2.5 mm

<u>TTA Tray</u>	For Forks, Cages & Plates - Empty	\$150.00
------------------------	-----------------------------------	----------

A tray will hold the following:

Plates

- 2 each of 3, 4, 5, 6, 7, & 8 hole

Forks

- 2 each of 3, 4, 5, 6, 7, & 8 prong


Cages

- 6 mm - 1 each of 16, 19, & 22 mm
- 9 mm - 1 each of 19, 22, & 25 mm
- 12 mm - 1 each of 22, 25, & 28 mm

<i>TTA Kit - Stainless</i>	Full kit with case, tray, and all items included, stainless implants
<i>TTA Kit - Titanium</i>	Full kit with case, tray, and all items included, titanium implants

Call for an estimate!

T.T.A. Kit



Code:

TTA Box B

Description:

Dimensions 280mm X 175mm X 95mm

Holds cages sizes 3, 4.5, 6, 7.5, 9, 10.5, 12, 15.
Holds plates and forks sizes 2, 3, 4, 5, 6, 7, 8.

Price for empty kit:

\$800.00

Contact us for pricing on a filled kit.



T.T.A. R Cages



New From V.O.I.

RTTA C3/10 T	TTA R Cage 3/10mm - Titanium	\$140.00
RTTA C3/13 T	TTA R Cage 3/13mm - Titanium	\$140.00
RTTA C3/16 T	TTA R Cage 3/16mm - Titanium	\$140.00
RTTA C4.5/12 T	TTA R Cage 4.5/12mm - Titanium	\$140.00
RTTA C4.5/15 T	TTA R Cage 4.5/15mm - Titanium	\$140.00
RTTA C4.5/18 T	TTA R Cage 4.5/18mm - Titanium	\$140.00
RTTA C6/16 T	TTA R Cage 6/16mm - Titanium	\$145.00
RTTA C6/19 T	TTA R Cage 6/19mm - Titanium	\$145.00
RTTA C6/22 T	TTA R Cage 5/22mm - Titanium	\$145.00
RTTA C7.5/16 T	TTA R Cage 7.5/16mm - Titanium	\$145.00
RTTA C7.5/19 T	TTA R Cage 7.5/19mm - Titanium	\$145.00
RTTA C7.5/22 T	TTA R Cage 7.5/22mm - Titanium	\$145.00
RTTA C9/19 T	TTA R Cage 9/19mm - Titanium	\$145.00
RTTA C9/22 T	TTA R Cage 9/22mm - Titanium	\$145.00
RTTA C9/25 T	TTA R Cage 9/25mm - Titanium	\$145.00
RTTA C10.5/19 T	TTA R Cage 10.5/19mm - Titanium	\$150.00
RTTA C10.5/22 T	TTA R Cage 10.5/22mm - Titanium	\$150.00
RTTA C10.5/25 T	TTA R Cage 10.5/25mm - Titanium	\$150.00
RTTA C12/22 T	TTA R Cage 12/22mm - Titanium	\$150.00
RTTA C12/25 T	TTA R Cage 12/25mm - Titanium	\$150.00
RTTA C12/28 T	TTA R Cage 12/28mm - Titanium	\$150.00



T.T.A. Implants - Stainless

Stainless Cages

Code:	Length:	Price:	Code:	Length:	Price:
TTA C3/10 S	10 mm	\$65.00	TTA C3/10 SL	10 mm	\$65.00
TTA C3/13 S	13 mm	\$65.00	TTA C3/10 SR	10 mm	\$65.00
TTA C3/16 S	16 mm	\$65.00	TTA C3/13 SL	13 mm	\$65.00
TTA C3/19 S	19 mm	\$65.00	TTA C3/13 SR	13 mm	\$65.00
			TTA C3/16 SL	16 mm	\$65.00
TTA C4.5/10 S	10 mm	\$65.00	TTA C3/16 SR	16 mm	\$65.00
TTA C4.5/12 S	12 mm	\$65.00			
TTA C4.5/13 S	13 mm	\$65.00	TTA C6/16 SL	16 mm	\$65.00
TTA C4.5/15 S	15 mm	\$65.00	TTA C6/16 SR	16 mm	\$65.00
TTA C4.5/16 S	16 mm	\$65.00	TTA C6/19 SL	19 mm	\$65.00
TTA C4.5/18 S	18 mm	\$65.00	TTA C6/19 SR	19 mm	\$65.00
			TTA C6/22 SL	22 mm	\$65.00
TTA C6/12 S	12 mm	\$65.00	TTA C6/22 SR	22 mm	\$65.00
TTA C6/13 S	13 mm	\$65.00			
TTA C6/16 S	16 mm	\$65.00	TTA C9/19 SL	19 mm	\$75.00
TTA C6/19 S	19 mm	\$65.00	TTA C9/19 SR	19 mm	\$75.00
TTA C6/22 S	22 mm	\$65.00	TTA C9/22 SL	22 mm	\$75.00
			TTA C9/22 SR	22 mm	\$75.00
TTA C7.5/16 S	16 mm	\$75.00	TTA C9/25 SL	25 mm	\$75.00
TTA C7.5/19 S	19 mm	\$75.00	TTA C9/25 SR	25 mm	\$75.00
TTA C7.5/22 S	22 mm	\$75.00			
			TTA C12/22 SL	22 mm	\$85.00
TTA C9/19 S	19 mm	\$75.00	TTA C12/22 SR	22 mm	\$85.00
TTA C9/22 S	22 mm	\$75.00	TTA C12/25 SL	25 mm	\$85.00
TTA C9/25 S	25 mm	\$75.00	TTA C12/25 SR	25 mm	\$85.00
			TTA C12/28 SL	28 mm	\$85.00
TTA C10.5/19 S	19 mm	\$85.00	TTA C12/28 SR	28 mm	\$85.00
TTA C10.5/22 S	22 mm	\$85.00			
TTA C10.5/25 S	25 mm	\$85.00	TTA C15/25 SL	25 mm	\$85.00
			TTA C15/25 SR	25 mm	\$85.00
TTA C12/22 S	22 mm	\$85.00	TTA C15/28 SL	28 mm	\$85.00
TTA C12/25 S	25 mm	\$85.00	TTA C15/28 SR	28 mm	\$85.00
TTA C12/28 S	28 mm	\$85.00	TTA C15/31 SL	31 mm	\$85.00
			TTA C15/31 SR	31 mm	\$85.00
TTA C15/25 S	25 mm	\$85.00			
TTA C15/28 S	28 mm	\$85.00			
TTA C15/31 S	31 mm	\$85.00			



Standard



Left



Right

“While using this product laterally, the distal edge of the cage is removed so that the anatomical fit is greater. There is no distal part of the cage sticking out and irritating the soft tissues surrounding the cage.”

— Radboud Kemme
Orthopedic Veterinarian
The Netherlands

T.T.A. Implants - Stainless

Stainless Cuttable Cages



Code:	Length:	Price:
TTA C6/22 S CUT	22mm	\$65.00
TTA C7.5/22 S CUT	22mm	\$75.00
TTA C9/26 S CUT	26mm	\$75.00
TTA C10.5/26 S CUT	26mm	\$85.00
TTA C12/28 S CUT	28mm	\$85.00
TTA C16/32 S CUT	32mm	\$95.00

Stainless Spacers



Code:	Width:	Thick:	Price:
S Spacer 2mm	10mm	2mm	\$15.60
S Spacer 4mm	10mm	4mm	\$15.60
S Spacer 6mm	10mm	6mm	\$15.60
S Spacer 9mm	10mm	9mm	\$15.60

T.T.A. Implants & Screws - Stainless

Stainless Forks

Code:	Prongs:	Price:
TTA F2 S	2	\$31.00
TTA F3 S	3	\$31.00
TTA F4 S	4	\$31.00
TTA F5 S	5	\$32.00
TTA F6 S	6	\$32.00
TTA F7 S	7	\$33.00
TTA F8 S	8	\$33.00



Code:	Length:	Holes:	Price:
TTA P2 S	41mm	2	\$27.00
TTA P3 S	42mm	3	\$27.00
TTA P4 S	52mm	4	\$27.00
TTA P5 S	65mm	5	\$27.00
TTA P6 S	78mm	6	\$35.00
TTA P7 S	91mm	7	\$35.00
TTA P8 S	104mm	8	\$35.00

Stainless Plates



TTA P3 S Long	48mm	3	\$27.00
TTA P4 S Long	58mm	4	\$27.00
TTA P5 S Long	71mm	5	\$27.00
TTA P6 S Long	84mm	6	\$35.00
TTA P7 S Long	97mm	7	\$35.00
TTA P8 S Long	110mm	8	\$35.00

Stainless Forkless Plates

Code:	Length	Holes:	Price:
TTA P2-2 S*	41mm	2	\$27.00
TTA P3-2 S*	42mm	2	\$27.00
TTA P4-2 S*	52mm	2	\$27.00
TTA P5-2 S*	65mm	2	\$27.00
TTA P6-3 S*	78mm	3	\$27.00
TTA P7-3 S*	91mm	3	\$35.00
TTA P8-3 S*	104mm	3	\$35.00



TTA P5-2 S Long*	71mm	2	\$27.00
TTA P6-3 S Long*	84mm	3	\$35.00
TTA P7-3 S Long*	97mm	3	\$35.00
TTA P8-3 S Long*	110mm	3	\$35.00



Patient Size (KG)

Less than 30 kg
30 - 45 kg
Over 45 kg

Patient Size (LB)

Less than 72 lbs
73 - 99 lbs
Over 100 lbs

Cage Size

6 mm
9 mm
12 mm

Plate Size

3 or 4
5 or 6
7 or 8

T.T.A. Screws - Stainless

2.4mm Cruciform Head Screws

Cruciform Head



Cruciform Head

Dia. of thread 2.4 mm
Thread pitch 1.0 mm
Dia. of core 1.7 mm
Dia. of head 4.0 mm

Drill for:

Threaded hole 1.8 mm
Gliding hole 2.4 mm

Code:	Length:	Price:
REG 240.06	6 mm	\$7.50
REG 240.08	8	\$7.50
REG 240.10	10	\$7.50
REG 240.12	12	\$7.75
REG 240.14	14	\$7.75
REG 240.16	16	\$7.75
REG 240.18	18	\$7.75
REG 240.20	20	\$7.75
REG 240.22	22	\$8.00
REG 240.24	24	\$8.00
REG 240.26	26	\$8.00
REG 240.28	28	\$8.00
REG 240.30	30	\$8.00
REG 240.32	32	\$8.25
REG 240.34	34	\$8.25
REG 240.36	36	\$8.25
REG 240.38	38	\$8.25
REG 240.40	40	\$8.25

2.4 mm Cruciform Head - Self-Tapping

Cruciform Head



Cruciform Head

Dia. of thread 2.4 mm
Thread pitch 1.0 mm
Dia. of core 1.7 mm
Dia. of head 4.0 mm

Drill for:

Threaded hole 1.8 mm
Gliding hole 2.4 mm

Code:	Length:	Price:
ST 240.06	6 mm	\$10.00
ST 240.08	8	\$10.00
ST 240.10	10	\$10.00
ST 240.12	12	\$10.00
ST 240.14	14	\$10.00
ST 240.16	16	\$10.00
ST 240.18	18	\$10.00
ST 240.20	20	\$10.00
ST 240.22	22	\$10.25
ST 240.24	24	\$10.25
ST 240.26	26	\$10.25
ST 240.28	28	\$10.25
ST 240.30	30	\$10.25
ST 240.32	32	\$10.50
ST 240.34	34	\$10.50
ST 240.36	36	\$10.50
ST 240.38	38	\$10.50
ST 240.40	40	\$10.50

2.4 mm Hex Head - Self-Tapping

Hex Head



Hex Head

Dia. of thread 2.4 mm
Thread pitch 1.0 mm
Dia. of core 1.7 mm
Dia. of head 4.0 mm

Drill for:

Threaded hole 1.8 mm
Gliding hole 2.4 mm

Code:	Length:	Price:
240.06T	6 mm	\$10.35
240.08T	8	\$10.35
240.10T	10	\$10.35
240.12T	12	\$10.64
240.14T	14	\$10.64
240.16T	16	\$10.64
240.18T	18	\$10.64
240.20T	20	\$10.64
240.22T	22	\$10.93
240.24T	24	\$10.93
240.26T	26	\$10.93
240.28T	28	\$10.93
240.30T	30	\$10.93
240.32T	32	\$11.21
240.34T	34	\$11.21
240.36T	36	\$11.21
240.38T	38	\$11.21
240.40T	40	\$11.21

T.T.A. Screws - Stainless

2.7 mm Hex Head - Self-Tapping

Code:	Length:	Price:
ST 270.10	10 mm	\$10.00
ST 270.12	12	\$10.25
ST 270.14	14	\$10.25
ST 270.16	16	\$10.25
ST 270.18	18	\$10.25
ST 270.20	20	\$10.25
ST 270.22	22	\$10.50
ST 270.24	24	\$10.50
ST 270.26	26	\$10.50
ST 270.28	28	\$10.50
ST 270.30	30	\$10.50
ST 270.32	32	\$10.75
ST 270.34	34	\$10.75

Hex Head



Hex Head

Dia. of thread	2.7 mm
Thread pitch	1.0 mm
Dia. of core	1.9 mm
Dia. of head	5.0 mm

Drill for:

Threaded hole	2.0 mm
Gliding hole	2.7 mm

For additional lengths, see page 2.4

3.5 mm Hex Head - Self-Tapping

Code:	Length:	Price:
ST 350.10	10 mm	\$9.75
ST 350.12	12	\$9.75
ST 350.14	14	\$9.75
ST 350.16	16	\$9.75
ST 350.18	18	\$9.75
ST 350.20	20	\$9.75
ST 350.22	22	\$10.00
ST 350.24	24	\$10.00
ST 350.26	26	\$10.00
ST 350.28	28	\$10.00
ST 350.30	30	\$10.00
ST 350.32	32	\$10.25
ST 350.34	34	\$10.25
ST 350.36	36	\$10.25
ST 350.38	38	\$10.25
ST 350.40	40	\$10.25

Hex Head



Hex Head

Dia. of thread	3.5 mm
Thread pitch	1.25 mm
Dia. of core	2.4 mm
Dia. of head	6.0 mm

Drill for:

Threaded hole	2.5 mm
Gliding hole	3.5 mm

For additional lengths, see page 2.5

T.T.A. Implants - Titanium

Titanium Cages

Code:	Length:	Price:	Code:	Length:	Price:
TTA C3/10 T	10 mm	\$75.00	TTA C3/10 TL	10 mm	\$75.00
TTA C3/13 T	13 mm	\$75.00	TTA C3/10 TR	10 mm	\$75.00
TTA C3/16 T	16 mm	\$75.00	TTA C3/13 TL	13 mm	\$75.00
TTA C3/19 T	19 mm	\$75.00	TTA C3/13 TR	13 mm	\$75.00
			TTA C3/16 TL	16 mm	\$75.00
TTA C4.5/10 T	10 mm	\$75.00	TTA C3/16 TR	16 mm	\$75.00
TTA C4.5/12 T	12 mm	\$75.00			
TTA C4.5/13 T	13 mm	\$75.00	TTA C6/16 TL	16 mm	\$75.00
TTA C4.5/15 T	15 mm	\$75.00	TTA C6/16 TR	16 mm	\$75.00
TTA C4.5/16 T	16 mm	\$75.00	TTA C6/19 TL	19 mm	\$75.00
TTA C4.5/18 T	18 mm	\$75.00	TTA C6/19 TR	19 mm	\$75.00
			TTA C6/22 TL	22 mm	\$75.00
TTA C6/13 T	13 mm	\$75.00	TTA C6/22 TR	22 mm	\$75.00
TTA C6/16 T	16 mm	\$75.00			
TTA C6/19 T	19 mm	\$75.00	TTA C9/19 TL	19 mm	\$85.00
			TTA C9/19 TR	19 mm	\$85.00
			TTA C9/22 TL	22 mm	\$85.00
TTA C7.5/16 T	16 mm	\$85.00	TTA C9/22 TR	22 mm	\$85.00
TTA C7.5/19 T	19 mm	\$85.00	TTA C9/25 TL	25 mm	\$85.00
			TTA C9/25 TR	25 mm	\$85.00
TTA C9/16 T	16 mm	\$85.00	TTA C12/22 TL	22 mm	\$95.00
TTA C9/19 T	19 mm	\$85.00	TTA C12/22 TR	22 mm	\$95.00
TTA C9/25 T	25 mm	\$85.00	TTA C12/25 TL	25 mm	\$95.00
			TTA C12/25 TR	25 mm	\$95.00
			TTA C12/28 TL	28 mm	\$95.00
			TTA C12/28 TR	28 mm	\$95.00
TTA C10.5/19 T	19 mm	\$95.00			
TTA C10.5/22 T	22 mm	\$95.00	TTA C15/25 TL	25 mm	\$95.00
TTA C10.5/25 T	25 mm	\$95.00	TTA C15/25 TR	25 mm	\$95.00
			TTA C15/28 TL	28 mm	\$95.00
TTA C12/22 T	22 mm	\$95.00	TTA C15/28 TR	28 mm	\$95.00
TTA C12/25 T	25 mm	\$95.00	TTA C15/31 TL	31 mm	\$95.00
TTA C12/28 T	28 mm	\$95.00	TTA C15/31 TR	31 mm	\$95.00
TTA C15/25 T	25 mm	\$95.00			
TTA C15/28 T	28 mm	\$95.00			
TTA C15/31 T	31 mm	\$95.00			



Standard



Left



Right

T.T.A. Implants - Titanium

Titanium Cuttable Cages

Code:	Length:	Price:
TTA C6/22 T CUT	22 mm	\$75.00
TTA C7.5/22 T CUT	22 mm	\$85.00
TTA C9/26 T CUT	26 mm	\$85.00
TTA C10.5/26 T CUT	26 mm	\$95.00
TTA C12/28 T CUT	28 mm	\$95.00
TTA C16/32 T CUT	32 mm	\$105.00



Titanium Spacers



Code:	Width:	Thick:	Price:
TI Spacer 2mm	10mm	2mm	\$15.60
TI Spacer 4mm	10mm	4mm	\$15.60
TI Spacer 6mm	10mm	6mm	\$15.60
TI Spacer 9mm	10mm	10mm	\$15.60

T.T.A. Implants - Titanium

Titanium Forks



Code:	Prongs:	Price:
TTA F2 T	2	\$33.00
TTA F3 T	3	\$33.00
TTA F4 T	4	\$33.00
TTA F5 T	5	\$36.00
TTA F6 T	6	\$36.00
TTA F7 T	7	\$40.00
TTA F8 T	8	\$40.00

Titanium Plates



Code:	Length:	Holes:	Price:
TTA P2 T	41 mm	2	\$30.00
TTA P3 T	42 mm	3	\$30.00
TTA P4 T	52 mm	4	\$33.00
TTA P5 T	65 mm	5	\$33.00
TTA P6 T	78 mm	6	\$36.00
TTA P7 T	91 mm	7	\$36.00
TTA P8 T	104 mm	8	\$36.00
TTA P3 T Long	48 mm	3	\$33.00
TTA P4 T Long	58 mm	4	\$33.00
TTA P5 T Long	71 mm	5	\$33.00
TTA P6 T Long	84 mm	6	\$36.00
TTA P7 T Long	97 mm	7	\$36.00
TTA P8 T Long	110 mm	8	\$36.00

Titanium Forkless Plates



Code:	Length:	Holes:	Price:
TTA P2-2 T*	41 mm	2	\$30.00
TTA P3-2 T*	42 mm	3	\$30.00
TTA P4-2 T*	52 mm	4	\$33.00
TTA P5-2 T*	65 mm	5	\$33.00
TTA P6-3 T*	78 mm	6	\$36.00
TTA P7-3 T*	91 mm	7	\$36.00
TTA P8-3 T*	104 mm	8	\$36.00
TTA P5-2 T Long*	71 mm	5	\$33.00
TTA P6-3 T Long*	84 mm	6	\$36.00
TTA P7-3 T Long*	97 mm	7	\$36.00
TTA P8-3 T Long*	110 mm	8	\$36.00



Patient Size (KG)

Less than 30 kg
30 - 45 kg
Over 45 kg

Patient Size (LB)

Less than 72 lbs
73 - 99 lbs
Over 100 lbs

Cage Size

6 mm
9 mm
12 mm

Plate Size

3 or 4
5 or 6
7 or 8

T.T.A. Screws - Titanium

Code:	Length:	Price:
TI REG 240.06	6mm	\$12.00
TI REG 240.08	8	\$12.00
TI REG 240.10	10	\$12.00
TI REG 240.12	12	\$12.00
TI REG 240.14	14	\$12.00
TI REG 240.16	16	\$12.00
TI REG 240.18	18	\$12.00
TI REG 240.20	20	\$12.00
TI REG 240.22	22	\$12.00
TI REG 240.24	24	\$12.00
TI REG 240.26	26	\$12.00
TI REG 240.28	28	\$12.00
TI REG 240.30	30	\$12.00
TI REG 240.32	32	\$12.00
TI REG 240.34	34	\$12.00
TI REG 240.36	36	\$12.00
TI REG 240.38	38	\$12.00
TI REG 240.40	40	\$12.00

Code:	Length:	Price:
TI REG 270.10	10mm	\$11.00
TI REG 270.12	12	\$11.00
TI REG 270.14	14	\$11.00
TI REG 270.16	16	\$11.00
TI REG 270.18	18	\$11.00
TI REG 270.20	20	\$11.00
TI REG 270.22	22	\$11.00
TI REG 270.24	24	\$11.00
TI REG 270.26	26	\$11.00
TI REG 270.28	28	\$11.00
TI REG 270.30	30	\$11.00
TI REG 270.32	32	\$11.00
TI REG 270.34	34	\$11.00

Code:	Length:	Price:
TI REG 350.10	10mm	\$10.00
TI REG 350.12	12	\$10.00
TI REG 350.14	14	\$10.00
TI REG 350.16	16	\$10.00
TI REG 350.18	18	\$10.00
TI REG 350.20	20	\$10.00
TI REG 350.22	22	\$10.00
TI REG 350.24	24	\$10.00
TI REG 350.26	26	\$10.00
TI REG 350.28	28	\$10.00
TI REG 350.30	30	\$10.00
TI REG 350.32	32	\$10.00
TI REG 350.34	34	\$10.00
TI REG 350.36	36	\$10.00
TI REG 350.38	38	\$10.00
TI REG 350.40	40	\$10.00

2.4 mm Cruciform Head

Cruciform Head



Cruciform Head

Dia. of thread	2.4 mm
Thread pitch	1.0 mm
Dia. of core	1.7 mm
Dia. of head	4.0 mm

Drill for:

Threaded hole	1.8 mm
Gliding hole	2.4 mm

Check for availability

2.7 mm Hex Head

Hex Head



Hex Head

Dia. of thread	2.7 mm
Thread pitch	1.0 mm
Dia. of core	1.9 mm
Dia. of head	5.0 mm

Drill for:

Threaded hole	2.0 mm
Gliding hole	2.7 mm

3.5 mm Hex Head

Hex Head



Hex Head

Dia. of thread	3.5 mm
Thread pitch	1.25 mm
Dia. of core	2.4 mm
Dia. of head	6.0 mm

Drill for:

Threaded hole	2.5 mm
Gliding hole	3.5 mm

T.T.A. Screws - Titanium

Self-Tapping

2.4 mm Cruciform Head

Cruciform Head



Cruciform Head

Dia. of thread 2.4 mm
 Thread pitch 1.0 mm
 Dia. of core 1.7 mm
 Dia. of head 4.0 mm

Drill for:

Threaded hole 1.8 mm
 Gliding hole 2.4 mm

Code:

TI ST 240.06
 TI ST 240.08
 TI ST 240.10
 TI ST 240.12
 TI ST 240.14
 TI ST 240.16
 TI ST 240.18
 TI ST 240.20
 TI ST 240.22
 TI ST 240.24
 TI ST 240.26
 TI ST 240.28
 TI ST 240.30
 TI ST 240.32
 TI ST 240.34
 TI ST 240.36
 TI ST 240.38
 TI ST 240.40

Length:

6mm
 8
 10
 12
 14
 16
 18
 20
 22
 24
 26
 28
 30
 32
 34
 36
 38
 40

Price:

\$18.00
 \$18.00
 \$18.00
 \$18.00
 \$18.00
 \$18.00
 \$18.00
 \$18.00
 \$18.00
 \$18.00
 \$18.00
 \$18.00
 \$18.00
 \$18.00
 \$18.00
 \$18.00
 \$18.00
 \$18.00

2.4 mm Hex Head

Hex Head



Hex Head

Dia. of thread 2.4 mm
 Thread pitch 1.0 mm
 Dia. of core 1.7 mm
 Dia. of head 4.0 mm

Drill for:

Threaded hole 1.8 mm
 Gliding hole 2.4 mm

Code:

TI 240.10T
 TI 240.12T
 TI 240.14T
 TI 240.16T
 TI 240.18T
 TI 240.20T
 TI 240.22T
 TI 240.24T
 TI 240.26T
 TI 240.28T
 TI 240.30T
 TI 240.32T
 TI 240.34T
 TI 240.36T
 TI 240.38T
 TI 240.40T

Length:

10mm
 12
 14
 16
 18
 20
 22
 24
 26
 28
 30
 32
 34
 36
 38
 40

Price:

\$15.00
 \$15.00
 \$15.00
 \$15.00
 \$15.00
 \$15.00
 \$15.00
 \$15.00
 \$15.00
 \$15.00
 \$15.00
 \$15.00
 \$15.00
 \$15.00
 \$15.00
 \$15.00
 \$15.00

T.T.A. Screws - Titanium

Self-Tapping

Code:	Length:	Price:
TI ST 270.10	10mm	\$16.00
TI ST 270.12	12	\$16.00
TI ST 270.14	14	\$16.00
TI ST 270.16	16	\$16.00
TI ST 270.18	18	\$16.00
TI ST 270.20	20	\$16.00
TI ST 270.22	22	\$16.00
TI ST 270.24	24	\$16.00
TI ST 270.26	26	\$16.00
TI ST 270.28	28	\$16.00
TI ST 270.30	30	\$16.00
TI ST 270.32	32	\$16.00
TI ST 270.34	34	\$16.00
TI ST 270.36	36	\$16.00
TI ST 270.38	38	\$16.00
TI ST 270.40	40	\$16.00

Code:	Length:	Price:
TI ST 350.10	10mm	\$15.00
TI ST 350.12	12	\$15.00
TI ST 350.14	14	\$15.00
TI ST 350.16	16	\$15.00
TI ST 350.18	18	\$15.00
TI ST 350.20	20	\$15.00
TI ST 350.22	22	\$15.00
TI ST 350.24	24	\$15.00
TI ST 350.26	26	\$15.00
TI ST 350.28	28	\$15.00
TI ST 350.30	30	\$15.00
TI ST 350.32	32	\$15.00
TI ST 350.34	34	\$15.00
TI ST 350.36	36	\$15.00
TI ST 350.38	38	\$15.00
TI ST 350.40	40	\$15.00

2.7 mm Hex Head

Hex Head



Hex Head

Dia. of thread	2.7 mm
Thread pitch	1.0 mm
Dia. of core	1.9 mm
Dia. of head	5.0 mm

Drill for:

Threaded hole	2.0 mm
Gliding hole	2.7 mm

3.5 mm Hex Head

Hex Head



Hex Head

Dia. of thread	3.5 mm
Thread pitch	1.25 mm
Dia. of core	2.4 mm
Dia. of head	6.0 mm

Drill for:

Threaded hole	2.5 mm
Gliding hole	3.5 mm

Cancellous

Code:	Length:	Price:
TI CAN 400.18	18mm	\$15.00
TI CAN 400.20	20	\$15.00
TI CAN 400.22	22	\$15.00
TI CAN 400.24	24	\$15.00
TI CAN 400.26	26	\$15.00

Ask for availability on additional lengths

4.0mm Cancellous Screw

Hex Head



Hex Head

Dia. of thread	4.0 mm
Thread pitch	1.75 mm
Dia. of core	1.9 mm
Dia. of head	6.0 mm

Drill for:

Threaded hole	2.5 mm
---------------	--------

T.T.A. Instruments



Code:	Description:	Price:
4.10.01	Drill Guide For The Forks, 8 Holes	\$145.00
4.10.01 A	Drill Guide For The Forks, 4 Holes	\$95.00



4.10.02	Fork Inserter	\$135.00
---------	---------------	----------



4.10.03	Plate Bender	\$125.00
---------	--------------	----------



4.10.04	T-Handle	\$125.00
---------	----------	----------



4.10.04 D	3.0 mm Spreader	\$10.00
-----------	-----------------	---------



4.10.04 H	4.5 mm Spreader	\$10.00
-----------	-----------------	---------



4.10.04 A	6.0 mm Spreader	\$10.00
-----------	-----------------	---------

4.10.04 F	7.5 mm Spreader	\$10.00
-----------	-----------------	---------

4.10.04 B	9.0 mm Spreader	\$10.00
-----------	-----------------	---------

4.10.04 G	10.5 mm Spreader	\$10.00
-----------	------------------	---------

4.10.04 C	12.0 mm Spreader	\$10.00
-----------	------------------	---------

4.10.04 E	15.0 mm Spreader	\$10.00
-----------	------------------	---------



4.10.05	1.9 mm Pins	\$7.00
---------	-------------	--------



4.20.01	Drill Sleeve 2.7/2.0	\$175.00
---------	----------------------	----------



4.20.02	Drill Sleeve 3.5 60 mm	\$175.00
---------	------------------------	----------



1.8 JAC 100	Drill Bit 1.8 mm	\$35.00
-------------	------------------	---------



2.0 JAC 85	Drill Bit 2.0 mm	\$35.00
------------	------------------	---------



2.5 JAC 95	Drill Bit 2.5 mm	\$35.00
------------	------------------	---------

T.T.A. Instruments

Code:	Description:	Price:
4.20.05	Depth Gauge - Fits 2.4 screws, measures up to 60mm	\$175.00
		
20020	Phenolic Handle with Quick Coupling	\$150.00
		
20021A	Qck Connect/AO Hex Shaft for 2.4 Screws	\$125.00
		
20022	Screwdriver Insert, 2.5 mm Hex Head	\$125.00
		
4.20.08	Screwdriver Insert, 2.4 mm Cruciform Head	\$125.00
		
10024	2.4mm Hex Screwdriver	\$161.00
		
20025 A	2.4mm Cruciform Screwdriver	\$161.00
		
20028	2.5 mm Hex Screwdriver for 2.7, 3.5, 4.0 screws	\$161.00
		

T.T.A. Instruments

Code: **Description:** **Price:**
 4.30.01 Hammer, 148g \$88.70



Code: **Description:** **Price:**
 46.2330 Bone Reduction Forceps, 5 1/4" \$130.38
 46.2350 Bone Reduction Forceps, 5" \$144.58
 46.2340 Bone Reduction Forceps, 8" \$167.06
 4.30.02 Bone Reduction Forceps, 5 1/4" \$104.10
 4.30.04 Bone Reduction Forceps, 8" \$137.00



Code: **Description:** **Price:**
 47.6190 Pelvic Reduction Forceps, 7 3/4" \$399.73
 47.6192 Pelvic Reduction Forceps, 9 1/2" \$408.95
 47.6196 Pelvic Reduction Forceps, 10" \$354.68



Code: **Description:** **Price:**
 Cage Cutter TTA Cage Cutter \$325.00

