

Titanium Screws

1.5 mm Cortex Head, Titanium

Code:	Length:	Price:
TI REG 150.06	6mm	\$10.00
TI REG 150.07	7	\$10.00
TI REG 150.08	8	\$10.00
TI REG 150.10	10	\$10.00
TI REG 150.11	11	\$10.00
TI REG 150.12	12	\$10.00
TI REG 150.13	13	\$10.00
TI REG 150.14	14	\$10.00
TI REG 150.15	15	\$10.00
TI REG 150.16	16	\$10.00
TI REG 150.17	17	\$10.00
TI REG 150.18	18	\$10.00
TI REG 150.19	19	\$10.00
TI REG 150.20	20	\$10.00



Cruciform Head

Cruciform Head

Dia. of thread	1.5 mm
Thread pitch	0.5 mm
Dia. of core	1.1 mm
Dia. of head	3.0 mm

Drill for:

Threaded hole	1.1 mm
Gliding hole	1.5 mm

Check for availability

2.0 mm Cortex Head, Titanium

Code:	Length:	Price:
TI REG 200.04	4mm	\$10.00
TI REG 200.05	5	\$10.00
TI REG 200.06	6	\$10.00
TI REG 200.07	7	\$10.00
TI REG 200.08	8	\$10.00
TI REG 200.09	9	\$10.00
TI REG 200.10	10	\$10.00
TI REG 200.11	11	\$10.00
TI REG 200.12	12	\$10.00
TI REG 200.13	13	\$10.00
TI REG 200.14	14	\$10.00
TI REG 200.15	15	\$10.00
TI REG 200.16	16	\$10.00
TI REG 200.17	17	\$10.00
TI REG 200.18	18	\$10.00
TI REG 200.19	19	\$10.00
TI REG 200.20	20	\$10.00



Cruciform Head

Cruciform Head

Dia. of thread	2.0 mm
Thread pitch	0.6 mm
Dia. of core	1.4 mm
Dia. of head	3.5 mm

Drill for:

Threaded hole	1.5 mm
Gliding hole	2.0 mm

Check for availability

2.4 mm Cortex Screws, Titanium

Code:	Length:	Price:
TI REG 240.06	6mm	\$12.00
TI REG 240.08	8	\$12.00
TI REG 240.10	10	\$12.00
TI REG 240.12	12	\$12.00
TI REG 240.14	14	\$12.00
TI REG 240.16	16	\$12.00
TI REG 240.18	18	\$12.00
TI REG 240.20	20	\$12.00
TI REG 240.22	22	\$12.00
TI REG 240.24	24	\$12.00
TI REG 240.26	26	\$12.00
TI REG 240.28	28	\$12.00
TI REG 240.30	30	\$12.00
TI REG 240.32	32	\$12.00
TI REG 240.34	34	\$12.00
TI REG 240.36	36	\$12.00
TI REG 240.38	38	\$12.00
TI REG 240.40	40	\$12.00



Cruciform Head

Cruciform Head

Dia. of thread	2.4 mm
Thread pitch	1.0 mm
Dia. of core	1.7 mm
Dia. of head	4.0 mm

Drill for:

Threaded hole	1.8 mm
Gliding hole	2.4 mm

Check for availability

Titanium Screws

2.7 mm Cortex Screws

Hex Head



Hex Head

Dia. of thread 2.7 mm
 Thread pitch 1.0 mm
 Dia. of core 1.9 mm
 Dia. of head 5.0 mm

Drill for:

Threaded hole 2.0 mm
 Gliding hole 2.7 mm

Code:

TI REG 270.10
 TI REG 270.12
 TI REG 270.14
 TI REG 270.16
 TI REG 270.18
 TI REG 270.20
 TI REG 270.22
 TI REG 270.24
 TI REG 270.26
 TI REG 270.28
 TI REG 270.30
 TI REG 270.32
 TI REG 270.34

Length:

10mm
 12
 14
 16
 18
 20
 22
 24
 26
 28
 30
 32
 34

Price:

\$11.00
 \$11.00
 \$11.00
 \$11.00
 \$11.00
 \$11.00
 \$11.00
 \$11.00
 \$11.00
 \$11.00
 \$11.00
 \$11.00
 \$11.00

3.5 mm Cortex Screws

Hex Head



Hex Head

Dia. of thread 3.5 mm
 Thread pitch 1.25 mm
 Dia. of core 2.4 mm
 Dia. of head 6.0 mm

Drill for:

Threaded hole 2.5 mm
 Gliding hole 3.5 mm

Code:

TI REG 350.10
 TI REG 350.12
 TI REG 350.14
 TI REG 350.16
 TI REG 350.18
 TI REG 350.20
 TI REG 350.22
 TI REG 350.24
 TI REG 350.26
 TI REG 350.28
 TI REG 350.30
 TI REG 350.32
 TI REG 350.34
 TI REG 350.36
 TI REG 350.38
 TI REG 350.40

Length:

10mm
 12
 14
 16
 18
 20
 22
 24
 26
 28
 30
 32
 34
 36
 38
 40

Price:

\$10.00
 \$10.00
 \$10.00
 \$10.00
 \$10.00
 \$10.00
 \$10.00
 \$10.00
 \$10.00
 \$10.00
 \$10.00
 \$10.00
 \$10.00
 \$10.00
 \$10.00

Titanium Screws

2.7 mm Cortex Screws, Titanium

Code:	Length:	Price:
T270.05	5mm	\$22.00
T270.06	6	\$22.00
T270.07	7	\$22.00
T270.08	8	\$22.00
T270.09	9	\$22.00
T270.10	10	\$22.00
T270.12	12	\$22.00
T270.14	14	\$22.00
T270.16	16	\$22.00
T270.18	18	\$22.00
T270.20	20	\$22.00
T270.22	22	\$22.00
T270.24	24	\$22.00



Hex Head



Hex Head

Dia. of thread	2.7 mm
Thread pitch	1.0 mm
Dia. of core	1.9 mm
Dia. of head	5.0 mm
Hex Recess Width	2.5 mm

Check for availability

4.0mm Cancellous Screw, Titanium

Code:	Length:	Price:
TI CAN 400.18	18mm	\$11.00
TI CAN 400.20	20	\$11.00
TI CAN 400.22	22	\$11.00
TI CAN 400.24	24	\$11.00
TI CAN 400.26	26	\$11.00



Hex Head



Hex Head

Dia. of thread	4.0 mm
Thread pitch	1.75 mm
Dia. of core	1.9 mm
Dia. of head	6.0 mm

Drill for:

Threaded hole	2.5 mm
---------------	--------

Titanium Screws

2.0 mm Cortex Screws - Self-Tapping, Titanium

Cruciform Head



Cruciform Head

Dia. of thread 2.0 mm
 Thread pitch 0.6 mm
 Dia. of core 1.4 mm
 Dia. of head 3.5 mm

Code:	Length:	Price:
TI ST 200.04	4mm	\$28.00
TI ST 200.05	5	\$28.00
TI ST 200.06	6	\$28.00
TI ST 200.07	7	\$28.00
TI ST 200.08	8	\$28.00
TI ST 200.09	9	\$28.00
TI ST 200.10	10	\$28.00
TI ST 200.11	11	\$28.00
TI ST 200.12	12	\$28.00
TI ST 200.13	13	\$28.00
TI ST 200.14	14	\$28.00
TI ST 200.15	15	\$28.00
TI ST 200.16	16	\$28.00
TI ST 200.17	17	\$28.00
TI ST 200.18	18	\$28.00
TI ST 200.19	19	\$28.00
TI ST 200.20	20	\$28.00

2.4 mm Cortex Screws - Self-Tapping, Titanium

Cruciform Head



Cruciform Head

Dia. of thread 2.4 mm
 Thread pitch 1.0 mm
 Dia. of core 1.7 mm
 Dia. of head 4.0 mm

Drill for:

Threaded hole 1.8 mm
 Gliding hole 2.4 mm

Check for availability

Code:	Length:	Price:
TI ST 240.06	6mm	\$18.00
TI ST 240.08	8	\$18.00
TI ST 240.10	10	\$18.00
TI ST 240.12	12	\$18.00
TI ST 240.14	14	\$18.00
TI ST 240.16	16	\$18.00
TI ST 240.18	18	\$18.00
TI ST 240.20	20	\$18.00
TI ST 240.22	22	\$18.00
TI ST 240.24	24	\$18.00
TI ST 240.26	26	\$18.00
TI ST 240.28	28	\$18.00
TI ST 240.30	30	\$18.00
TI ST 240.32	32	\$18.00
TI ST 240.34	34	\$18.00
TI ST 240.36	36	\$18.00
TI ST 240.38	38	\$18.00
TI ST 240.40	40	\$18.00

Titanium Screws

2.7 mm Cortex Screw - Self-Tapping, Titanium

Code:	Length:	Price:
TI ST 270.10	10mm	\$16.00
TI ST 270.12	12	\$16.00
TI ST 270.14	14	\$16.00
TI ST 270.16	16	\$16.00
TI ST 270.18	18	\$16.00
TI ST 270.20	20	\$16.00
TI ST 270.22	22	\$16.00
TI ST 270.24	24	\$16.00
TI ST 270.26	26	\$16.00
TI ST 270.28	28	\$16.00
TI ST 270.30	30	\$16.00
TI ST 270.32	32	\$16.00
TI ST 270.34	34	\$16.00



Hex Head
 Dia. of thread 2.7 mm
 Thread pitch 1.0 mm
 Dia. of core 1.9 mm
 Dia. of head 5.0 mm

Drill for:
 Threaded hole 2.0 mm
 Gliding hole 2.7 mm

3.5 mm Cortex Screw - Self-Tapping, Titanium

Code:	Length:	Price:
TI ST 350.10	10mm	\$15.00
TI ST 350.12	12	\$15.00
TI ST 350.14	14	\$15.00
TI ST 350.16	16	\$15.00
TI ST 350.18	18	\$15.00
TI ST 350.20	20	\$15.00
TI ST 350.22	22	\$15.00
TI ST 350.24	24	\$15.00
TI ST 350.26	26	\$15.00
TI ST 350.28	28	\$15.00
TI ST 350.30	30	\$15.00
TI ST 350.32	32	\$15.00
TI ST 350.34	34	\$15.00
TI ST 350.36	36	\$15.00
TI ST 350.38	38	\$15.00
TI ST 350.40	40	\$15.00



Hex Head
 Dia. of thread 3.5 mm
 Thread pitch 1.25 mm
 Dia. of core 2.4 mm
 Dia. of head 6.0 mm

Drill for:
 Threaded hole 2.5 mm
 Gliding hole 3.5 mm

Titanium Screws

2.0 mm Cortex Screws - Self-Tapping, Titanium

Hex Head



Hex Head

Dia. of thread	2.0 mm
Thread pitch	0.6 mm
Dia. of core	1.4 mm
Dia. of head	4.0 mm
Hex Recess Width	1.5 mm

Code:

T200.05T
T200.06T
T200.07T
T200.08T
T200.10T
T200.12T
T200.14T

Length:

5mm
6
7
8
10
12
14

Price:

\$25.00
\$25.00
\$25.00
\$25.00
\$25.00
\$25.00
\$25.00

2.7 mm Cortex Screws - Self-Tapping, Titanium

Hex Head



Hex Head

Dia. of thread	2.7 mm
Thread pitch	1.0 mm
Dia. of core	1.9 mm
Dia. of head	5.0 mm
Hex Recess Width	2.5 mm

Code:

T270.05T
T270.06T
T270.07T
T270.08T
T270.10T
T270.12T
T270.14T
T270.16T
T270.18T
T270.20T
T270.24T

Length:

5mm
6
7
8
10
12
14
16
18
20
24

Price:

\$25.00
\$25.00
\$25.00
\$25.00
\$25.00
\$25.00
\$25.00
\$25.00
\$25.00
\$25.00

Check for availability

Titanium Screws

1.5 mm Star Head, Cortical Titanium - Self-Tapping, RU-6

Code:	Length:	Price:
05.16.06	6mm	\$19.84
05.16.07	7	\$19.84
05.16.08	8	\$19.84
05.16.10	10	\$19.84
05.16.12	12	\$19.84
05.16.14	14	\$19.84
05.16.16	16	\$19.84
05.16.18	18	\$19.84
05.16.20	20	\$19.84
05.16.22	22	\$19.84
05.16.24	24	\$19.84
05.16.26	26	\$19.84
05.16.28	28	\$19.84
05.16.30	30	\$19.84



Star Head



2.4 mm Locking Star Head, Titanium - Self-Tapping, RU-6

Star Head



Code:	Length:	Price:
05.15.06	6mm	\$19.84
05.15.08	8	\$19.84
05.15.10	10	\$19.84
05.15.12	12	\$19.84
05.15.14	14	\$19.84
05.15.16	16	\$19.84

2.4 mm Star Head, Cortical Titanium - Self-Tapping, RU-10

Code:	Length:	Price:
05.26.10	10mm	\$19.84
05.26.12	12	\$19.84
05.26.14	14	\$19.84
05.26.16	16	\$19.84
05.26.18	18	\$19.84
05.26.20	20	\$19.84
05.26.22	22	\$19.84
05.26.24	24	\$19.84
05.26.26	26	\$19.84
05.26.28	28	\$19.84
05.26.30	30	\$19.84
05.26.32	32	\$19.84



Star Head



Titanium Screws

3.2 mm Locking Star Head, Titanium - Self-Tapping, RU-10

Star Head



Code:	Length:	Price:
05.25.08	8mm	\$19.84
05.25.10	10	\$19.84
05.25.12	12	\$19.84
05.25.14	14	\$19.84
05.25.16	16	\$19.84
05.25.18	18	\$19.84

2.7 mm Star Head, Cortical Titanium - Self-Tapping, RU-10

Code:	Length:	Price:
05.36.12	12mm	\$19.84
05.36.14	14	\$19.84
05.36.16	16	\$19.84
05.36.18	18	\$19.84
05.36.20	20	\$19.84
05.36.22	22	\$19.84
05.36.24	24	\$19.84
05.36.26	26	\$19.84
05.36.28	28	\$19.84
05.36.30	30	\$19.84
05.36.32	32	\$19.84
05.36.34	34	\$19.84



Star Head



Titanium Screws

4.0 mm Locking Star Head, Titanium - Self-Tapping, RU-10

Code:	Length:	Price:
05.35.10	10mm	\$19.84
05.35.12	12	\$19.84
05.35.14	14	\$19.84
05.35.16	16	\$19.84
05.35.18	18	\$19.84
05.35.20	20	\$19.84
05.35.22	22	\$19.84
05.35.24	24	\$19.84



Star Head



3.5 mm Locking Star Head Screws - Self-Tapping

Code:	Length:	Price:
T SLH ST 350.10	10mm	\$25.00
T SLH ST 350.12	12	\$25.00
T SLH ST 350.14	14	\$25.00
T SLH ST 350.16	16	\$25.00
T SLH ST 350.18	18	\$25.00
T SLH ST 350.20	20	\$25.00
T SLH ST 350.22	22	\$25.00
T SLH ST 350.24	24	\$25.00
T SLH ST 350.26	26	\$25.00
T SLH ST 350.28	28	\$25.00
T SLH ST 350.30	30	\$25.00
T SLH ST 350.32	32	\$25.00
T SLH ST 350.34	34	\$25.00
T SLH ST 350.36	36	\$25.00
T SLH ST 350.38	38	\$25.00
T SLH ST 350.40	40	\$25.00
T SLH ST 350.42	42	\$25.00
T SLH ST 350.44	44	\$25.00
T SLH ST 350.45	45	\$25.00
T SLH ST 350.46	46	\$25.00
T SLH ST 350.48	48	\$25.00
T SLH ST 350.50	50	\$25.00
T SLH ST 350.52	52	\$25.00
T SLH ST 350.54	54	\$25.00
T SLH ST 350.56	56	\$25.00
T SLH ST 350.58	58	\$25.00
T SLH ST 350.60	60	\$25.00

Star Head

Dia of Thread 3.5mm
Thread pitch 0.8mm
Dia of core 2.9mm
Dia of head 5.0mm



Star Head



Drill for:
Threaded hole 2.8mm

Drill Guide for Locking Screws

Code	Price
LDG 3.5	\$125.00

Titanium Implants

Soft Tissue Bone Anchors - Titanium

Code:	Length:	Price:
TI AN 8	8mm	\$35.00
TI AN 13	13	\$35.00




Titanium Suture Anchor Posts

Code:	Dia.	Length:	Price:
E5608	3.5mm	18mm	\$45.00
E5609	3.5	20	\$45.00
E5610	3.5	22	\$45.00
E5611	3.5	24	\$45.00
E5660	2.7	12	\$42.00
E5661	2.7	14	\$42.00
E5662	2.7	16	\$42.00
E5663	2.7	18	\$42.00




1.5 mm Plates - Titanium

Code:	Description:	Holes:	Length:	Price:
TI 1.520	Cutttable	20	100mm	\$154.00
TI 1.5 3T9	T-Plate	9	50	\$133.00
TI 1.5 4T9	T-Plate	9	50	\$145.00



2.0 mm Plates - Titanium

Code:	Description:	Holes:	Length:	Price:
TI 2.020	Cutttable	20	100mm	\$153.00
TI 2.0 3T9	T-Plate	9	50	\$133.00
TI 2.0 4T9	T-Plate	9	50	\$145.00



2.7 mm Plates - Titanium

Code:	Description:	Holes:	Length:	Price:
TI 2.706 CP	Compression	6	52mm	\$193.75
TI 2.708 CP	Compression	8	68	\$206.25



Titanium Implants


O Style Non-locking T.P.L.O. Plates - Titanium 3.5 mm

Code:		Description:	Price:
TI O 3.5 DLC L	Left	Titanium Low Contact Compression Holes	95.00
TI O 3.5 DLC R	Right		



Y Style T.P.L.O. Plates - Titanium Medium

TI YML LC	Left	Titanium Low Contact	135.00
TI YMR LC	Right		



Shank Mini Interface Half Pin, Titanium

Code:	Length:	Price:
I03035 T	0.035"	\$18.00
I03062 T	0.062"	\$18.00

Titanium Implants

Pearl Plates - Titanium

Code:	Length:	Description:	Price:	
T Pearl * 2.012	98 mm	Titanium 2.0 mm 12 hole	\$108.00	
T Pearl 2.416	130 mm	Titanium 2.4 mm 16 hole	\$134.40	* compatible with PBI 2.4
T Pearl 2.716	130 mm	Titanium 2.7 mm 16 hole	\$192.00	* compatible with PBI 2.7
T Pearl 3.520	244 mm	Titanium 3.5 mm 20 hole	\$234.00	* compatible with PBI 3.5

Pearl Plugs - Titanium

Code	Description	Price
T *2.0 Plug	Titanium 2.0 Pearl Plug	\$30.00
T 2.4 Plug	Titanium 2.4 Pearl Plug	\$15.60
T 2.7 Plug	Titanium 2.7 Pearl Plug	\$15.60
T 3.5 Plug	Titanium 3.5 Pearl Plug	\$14.40