

V.O.I. Suture

V.O.I. distributes wound closure products from leading U.S. companies that produce suture for both the veterinary and human market. We pride our sutures with a 100% satisfaction guarantee. If you're not satisfied with the suture for any reason, we will provide a credit, replace the product, or issue your money back.

Ethicon Trademark Name:

PDS®
VICRYL®
MONOCRYL®
PROLENE®


V.O.I. Material Equivalent

PDO
PGA
PGCL
Polypropylene

Needle Types and Code References


PRECISION REVERSE CUTTING NEEDLES

3/8 Circle




REVERSE CUTTING NEEDLES

3/8 Circle




TAPER POINT NEEDLES

1/2 Circle



1/2 Circle



STRAIGHT CUTTING NEEDLES

NKS 60mm



All Suture is sold 12 to a box

CP Medical	DemeTECH	Sutumed	Ethicon Needle Code	Material	Suture Length (Inches)
311A	PGA283022B0	ST PGA3022MR	SH-1	PGA	30"
339A	PGA282037B0	ST PGA3022MR	CT-1	PGA	30"
340A	PGA280037B0	ST PGA0037MR	CT-1	PGA	30"
385A	PGA185013F4	ST PGA5013TC18	C-3	PGA	18"
391A	PGA185019F4	ST PGA5019TC18	FS-2	PGA	18"
397A	PGA284019F4	ST PGA4019TC	FS-2	PGA	30"
398A	PGA283019F4	ST PGA3019TC	FS-2	PGA	30"
452A	PBA2863024F4	ST PGA3024TC	FS-1	PGA	30"
453A	PGA282024F4	ST PGA2024TC	FS-1	PGA	30"
467A	PGA280036B4	ST PGA0036MC	CP-1	PGA	30"
468A	PGA281036B4	ST PGA0136MC	CP-1	PGA	30"
492A	PGA186013F13	ST PGA6013TC18	P-3	PGA	18"
494A	PGA184013F13	ST PGA4013TC18	P-3	PGA	18"
586A	PGA282030F13	ST PGA2030TC	FSL	PGA	30"
587A	PGA280030F13	ST PGA0030TC	FSL	PGA	30"
637CG	CC232024F4	ST637CG	FS-1	Chromic Gut	30"
628B	NL262060A6	ST628B	KS	Nylon	30"
662B	NL164019F4	ST662B	FS-2	Nylon	18"
662M	NL294019F4	ST662M	FS-2	Nylon	30"
663B	NL163024F4	ST663B	FS-1	Nylon	18"
663M	NL293024F4	ST663M	FS-1	Nylon	30"
664B	NL162026F4	ST664B	FS	Nylon	18"
664M	NL292026F4	ST664M	FS	Nylon	30"
665M	NL292026F4		FS	Nylon	30"
1669B	NL163019F13		PS-2	Nylon	18"
928B	NL163019F4	ST928B	FS-2	Nylon	18"
684S	NL163024F4		FS-1	Silk	18"
685S	SK162026F4		FS	Silk	18"
632S	SK163023B4		X-A	Silk	18"
8411P	PM292026B0	ST8411P	CT-2	Polypropylene	30"
8423P	PM292037B0	ST8423P	CT-1	Polypropylene	30"
8424P	PM290037B0	ST8424P	CT-1	Polypropylene	30"
8683P	PM194019F4	ST8683P	FS-2	Polypropylene	18"
8684P	PM193024F4	ST8684P	FS-1	Polypropylene	18"
8685P	PM192026F4	ST8685P	FS	Polypropylene	18"
8686P	PM195019F13		PS-2	Polypropylene	18"
8689P	PM292040F4	ST8689P	FSLX	Polypropylene	30"
8832P	PM293026B0	ST8832P	SH	Polypropylene	30"
L310	MO284022B0		SH-1	PGCL	30"
L315	MO284026B0		SH	PGCL	30"
L316	MO283026B0		SH	PGCL	30"
L317	MO282026B0		SH	PGCL	30"

Color	Thread Dia.	Needle Size	Needle Curvature	Needle Style	Pack Type
Violet	3/0	22mm	1/2 Circle	Taper Point	12ct Box
Violet	2/0	37mm	1/2 Circle	Taper Point	12ct Box
Violet	0/0	37mm	1/2 Circle	Taper Point	12ct Box
Violet	5/0	13mm	3/8 Circle	Reverse Cutting	12ct Box
Violet	5/0	19mm	3/8 Circle	Reverse Cutting	12ct Box
Violet	4/0	19mm	3/8 Circle	Reverse Cutting	12ct Box
Violet	3/0	19mm	3/8 Circle	Reverse Cutting	12ct Box
Violet	3/0	24mm	3/8 Circle	Reverse Cutting	12ct Box
Violet	2/0	24mm	3/8 Circle	Reverse Cutting	12ct Box
Violet	0/0	36mm	1/2 Circle	Reverse Cutting	12ct Box
Violet	1/0	36mm	1/2 Circle	Reverse Cutting	12ct Box
Violet	6/0	13mm	3/8 Circle	Precision Reverse Cutting	12ct Box
Violet	4/0	13mm	3/8 Circle	Precision Reverse Cutting	12ct Box
Violet	2/0	30mm	3/8 Circle	Precision Reverse Cutting	12ct Box
Violet	0/0	30mm	3/8 Circle	Precision Reverse Cutting	12ct Box
Beige	2/0	24mm	3/8 Circle	Reverse Cutting	12ct Box
Black	2/0	60mm	Straight	Reverse Cutting	12ct Box
Black	4/0	19mm	3/8 Circle	Reverse Cutting	12ct Box
Blue	4/0	19mm	3/8 Circle	Reverse Cutting	12ct Box
Black	3/0	24mm	3/8 Circle	Reverse Cutting	12ct Box
Blue	3/0	24mm	3/8 Circle	Reverse Cutting	12ct Box
Black	2/0	26mm	3/8 Circle	Reverse Cutting	12ct Box
Blue	2/0	26mm	3/8 Circle	Reverse Cutting	12ct Box
Blue	2/0	26mm	3/8 Circle	Reverse Cutting	12ct Box
Black	3/0	19mm	3/8 Circle	Precision Reverse Cutting	12ct Box
Black	3/0	19mm	3/8 Circle	Reverse Cutting	12ct Box
Black	3/0	24mm	3/8 Circle	Reverse Cutting	12ct Box
Black	2/0	26mm	3/8 Circle	Reverse Cutting	12ct Box
Black	3/0	23mm	1/2 Circle	Reverse Cutting	12ct Box
Blue	2/0	26mm	1/2 Circle	Taper Point	12ct Box
Blue	2/0	37mm	1/2 Circle	Taper Point	12ct Box
Blue	0/0	37mm	1/2 Circle	Taper Point	12ct Box
Blue	4/0	19mm	3/8 Circle	Reverse Cutting	12ct Box
Blue	3/0	24mm	3/8 Circle	Reverse Cutting	12ct Box
Blue	2/0	26mm	3/8 Circle	Reverse Cutting	12ct Box
Blue	5/0	19mm	3/8 Circle	Precision Reverse Cutting	12ct Box
Blue	2/0	40mm	3/8 Circle	Reverse Cutting	12ct Box **
Blue	3/0	26mm	1/2 Circle	Taper Point	12ct Box
Violet	4/0	22mm	1/2 Circle	Taper Point	12ct Box
Violet	4/0	26mm	1/2 Circle	Taper Point	12ct Box
Violet	3/0	26mm	1/2 Circle	Taper Point	12ct Box
Violet	2/0	26mm	1/2 Circle	Taper Point	12ct Box **

** Built To Order CP Medical Only

CP Medical	DemeTECH	Sutumed	Ethicon Needle Code	Material	Suture Length (Inches)
L344	MO1083037B0		CT-1	PGCL	36"
L345	MO1082037B0		CT-1	PGCL	36"
L463	MO1080037B0		CT-1	PGCL	36"
L493	MO185013F13		P-3	PGCL	18"
L494	MO174013F13		P-3	PGCL	18"
L762	MO1082026B0		CT-2	PGCL	36"
L922	MO1084019F4		FS-2	PGCL	36"
L923	MO1083019F4		FS-2	PGCL	36"
L942	MO1083024F4		FS-1	PGCL	36"
L943	MO1082024F4		FS-1	PGCL	36"
L967	MO1080036B4		CP-1	PGCL	36"
L987	MO280030F4		FSL	PGCL	30"
M304	PX284017B0	STM304	RB-1	PDO	30"
M305	PX283017B0	STM305	RB-1	PDO	30"
M310	PX284022B0	STM310	SH-1	PDO	30"
M311	PX283022B0	STM311	SH-1	PDO	30"
M315	PX284026B0	STM315	SH	PDO	30"
M316	PX283026B0	STM316	SH	PDO	30"
M317	PX282026B0	STM317	SH	PDO	30"
M332	PX283026B0	STM332	CT-2	PDO	30"
M333	PX282026B0	STM333	CT-2	PDO	30"
M334	PX280026B0	STM334	CT-2	PDO	30"
M338	PX283037B0	STM338	CT-1	PDO	30"
M339	PX282037B0	STM339	CT-1	PDO	30"
M340	PX280037B0	STM340	CT-1	PDO	30"
M341	PX281037B0	STM341	CT-1	PDO	30"
M397	PX284019F4	STM397	FS-2	PDO	30"
M398	PX283019F4	STM398	FS-2	PDO	30"
M451	PX282024F4	STM451	FS-1	PDO	30"
M452	PX283024F4	STM452	FS-1	PDO	30"
M466	PX282036B4	STM466	CP-1	PDO	30"
M467	PX280036B4	STM467	CP-1	PDO	30"
M468	PX281036B4	STM468	CP-1	PDO	30"
M968	PX283026B0		CP-2	PDO	30"
M969	PX282026B0		CP-2	PDO	30"
M970	PX280026B0		CP-2	PDO	30"
8685GP	N/A		FS	GLO-STITCH™	30"
8684GP	N/A		FS-1	GLO-STITCH™	30"
8684FF	N/A		FS-1	FLUOROFIL®	30"
8685FF	N/A		FS	FLUOROFIL®	30"

Color	Tread Dia.	Needle Size	Needle Curvature	Needle Style	Pack Type
Violet	3/0	37mm	1/2 Circle	Taper Point	12ct Box
Violet	2/0	37mm	1/2 Circle	Taper Point	12ct Box
Violet	0/0	37mm	1/2 Circle	Taper Point	12ct Box
Violet	5/0	13mm	3/8 Circle	Precision Reverse Cutting	12ct Box
Undyed	4/0	13mm	3/8 Circle	Precision Reverse Cutting	12ct Box
Violet	2/0	26mm	1/2 Circle	Taper Point	12ct Box **
Violet	4/0	19mm	3/8 Circle	Reverse Cutting	12ct Box
Violet	3/0	19mm	3/8 Circle	Reverse Cutting	12ct Box
Violet	3/0	24mm	3/8 Circle	Reverse Cutting	12ct Box
Violet	2/0	24mm	3/8 Circle	Reverse Cutting	12ct Box
Violet	0/0	36mm	1/2 Circle	Reverse Cutting	12ct Box
Violet	0/0	30mm	3/8 Circle	Taper Point	12ct Box
Violet	4/0	17mm	1/2 Circle	Taper Point	12ct Box
Violet	3/0	17mm	1/2 Circle	Taper Point	12ct Box
Violet	4/0	22mm	1/2 Circle	Taper Point	12ct Box
Violet	3/0	22mm	1/2 Circle	Taper Point	12ct Box
Violet	4/0	26mm	1/2 Circle	Taper Point	12ct Box
Violet	3/0	26mm	1/2 Circle	Taper Point	12ct Box
Violet	2/0	26mm	1/2 Circle	Taper Point	12ct Box
Violet	3/0	26mm	1/2 Circle	Taper Point	12ct Box
Violet	2/0	26mm	1/2 Circle	Taper Point	12ct Box
Violet	0/0	26mm	1/2 Circle	Taper Point	12ct Box
Violet	3/0	37mm	1/2 Circle	Taper Point	12ct Box
Violet	2/0	37mm	1/2 Circle	Taper Point	12ct Box
Violet	0/0	37mm	1/2 Circle	Taper Point	12ct Box
Violet	1/0	37mm	1/2 Circle	Reverse Cutting	12ct Box
Violet	4/0	19mm	3/8 Circle	Reverse Cutting	12ct Box
Violet	3/0	19mm	3/8 Circle	Reverse Cutting	12ct Box
Violet	2/0	24mm	3/8 Circle	Reverse Cutting	12ct Box
Violet	3/0	24mm	3/8 Circle	Reverse Cutting	12ct Box
Violet	2/0	36mm	1/2 Circle	Reverse Cutting	12ct Box
Violet	0/0	36mm	1/2 Circle	Reverse Cutting	12ct Box
Violet	1/0	36mm	1/2 Circle	Taper Point	12ct Box
Violet	3/0	26mm	1/2 Circle	Taper Point	12ct Box
Violet	2/0	26mm	1/2 Circle	Taper Point	12ct Box
Violet	0/0	26mm	1/2 Circle	Reverse Cutting	12ct Box
Green	2/0	26mm	3/8 Circle	Reverse Cutting	12ct Box
Green	3/0	24mm	3/8 Circle	Reverse Cutting	12ct Box
Pink	3/0	24mm	3/8 Circle	Reverse Cutting	12ct Box
Pink	2/0	26mm	3/8 Circle	Reverse Cutting	12ct Box

** Built To Order CP Medical Only

Other sizes available.
Contact your V.O.I. account representative
today for more details.

V.O.I. Suture - Staple and Glue -

Stapler

CP35R	35 Regular	5.4mm x 3.6mm	\$12.00
CP35W	35 Wide	7.0mm x 4.0mm	\$12.00



Glue

CA1999	2ml	1 bottle, 12 applicator tips	\$15.00
--------	-----	------------------------------	---------

