

2.4 & 3.0 mm Countersink Compression Screws

Cancellous Partial Thread

2.4 mm Countersink Compression Screws - T8

Code:	Length:	Price:
CCS 240.18	18mm	\$70.00
CCS 240.20	20	\$70.00
CCS 240.22	22	\$70.00
CCS 240.24	24	\$70.00
CCS 240.26	26	\$70.00

Star Head



Diameter of thread:	2.0 mm
Thread pitch:	1.0 mm
Diameter of core:	2.0 mm
Diameter of head:	3.0 mm



3.0 mm Countersink Compression Screws - T8

Code:	Length:	Price:
CCS 300.20	20mm	\$70.00
CCS 300.22	22	\$70.00
CCS 300.24	24	\$70.00
CCS 300.26	26	\$70.00
CCS 300.28	28	\$70.00
CCS 300.30	30	\$70.00
CCS 300.32	32	\$70.00
CCS 300.34	34	\$70.00
CCS 300.36	36	\$70.00
CCS 300.38	38	\$70.00
CCS 300.40	40	\$70.00

Star Head



Diameter of thread:	3.0 mm
Thread pitch:	1.25 mm
Diameter of core:	2.0 mm
Diameter of head:	3.5 mm



4.5 mm Countersink Compression Screws

Cancellous Partial Thread

4.5 mm Countersink Compression Screws - T15

Code:	Length:	Price:
CCS 450.28	28mm	\$70.00
CCS 450.30	30	\$70.00
CCS 450.32	32	\$70.00
CCS 450.34	34	\$70.00
CCS 450.36	36	\$70.00
CCS 450.38	38	\$70.00
CCS 450.40	40	\$70.00
CCS 450.42	42	\$70.00
CCS 450.44	44	\$70.00
CCS 450.46	46	\$70.00
CCS 450.48	48	\$70.00
CCS 450.50	50	\$70.00
CCS 450.52	52	\$70.00
CCS 450.54	54	\$70.00
CCS 450.56	56	\$70.00
CCS 450.58	58	\$70.00
CCS 450.60	60	\$70.00
CCS 450.65	65	\$70.00
CCS 450.70	70	\$70.00
CCS 450.75	75	\$70.00

Star Head





Diameter of thread: 4.5 mm
Thread pitch: 1.75 mm
Diameter of core: 3.0 mm
Diameter of head: 5.0 mm



C.B.L.O.

CBLO Instrumentation

	Code:	Description:	Price:
	CCS-CS-2.4	Compression sleeve for 2.4mm compression screw	\$225.00
	CCS-CS-3.0	Compression sleeve for 3.0mm compression screw	\$225.00
	CCS-CS-4.5	Compression sleeve for 4.5mm compression screw	\$225.00

	Code:	Description:	Price:
	CCS-HDL-2.4/3.0	Compression sleeve handle for 2.4/3.0mm countersink compression screw	\$225.00
	CCS-HDL-4.5	Compression sleeve handle for 4.5mm countersink compression screw	\$225.00



Code:	Description:	Price:
CCS-CD-T8	Cannulated star drive screwdriver, T8, For CCS 2.4 & 3.0	\$267.00
CCS-CD-T15	Cannulated star drive screwdriver, T15, For CCS 4.5	\$267.00

Code:	Length:	Use With:	Price:
CAN 2.0 150 QCK	150mm	0.045 K Wire	\$198.00
CAN 3.2 170 QCK	170mm	0.054 K Wire	\$198.00
CAN 3.2 225 QCK	225mm	0.062 K Wire	\$198.00
2.5 JAC 130 CAN	130mm	0.032 K Wire	\$178.00
2.7 JAC 130 CAN	130mm	0.045 K Wire	\$178.00

Cannulated drill bits are used with this system.

C.B.L.O.

CBLO Plates

Code: DT L CBLO 2.0 L DT L CBLO 2.0 R	Features: Compression holes, locking, low contact	Use With: 0.027 K Wire	Price: \$120.00 \$120.00
Code: DT L CBLO 2.4 L DT L CBLO 2.4 R	Features: Compression holes, locking, low contact	0.027 K Wire	Price: \$120.00 \$120.00
Code: DT L CBLO 2.7 L DT L CBLO 2.7 R	Features: Compression holes, locking, low contact	0.045 K Wire	Price: \$120.00 \$120.00
Code: DT L CBLO 3.5 SM L DT L CBLO 3.5 SM R	Features: Compression holes, locking, low contact	0.045 K Wire	Price: \$125.00 \$125.00
Code: DT L CBLO 3.5 L DT L CBLO 3.5 R	Features: Compression holes, locking, low contact	0.045 K Wire	Price: \$125.00 \$125.00
Code: DT L CBLO 3.5B L DT L CBLO 3.5B R	Features: Compression holes, locking, low contact	0.045 K Wire	Price: \$135.00 \$135.00



DT L CBLO 3.5 Right and Left